Chemical Curing Adhesive

EN	Chemical curing adhesive
LT	Cheminio kietėjimo surišėjas
DE	Chemisch härtender Klebstoff
BG	Лепило за химическо втвърдяване
CS	Chemicky tuhnutí lepidlo
DA	Kemisk hærdende klæbemiddel
ET	Keemiliselt kõvenev liim
FR	Adhésif à durcissement chimique
EL	Κόλλα χημικής σκληρύνσης
HU	Kémiai kötésű ragasztó
IT	Adesivo a polimerizzazione chimica
LV	Ķīmiski cietējoša lieme
NO	Kjemisk herdende lim
PL	Klej utwardzany chemicznie
PT	Adesivo de cura química
RO	Adeziv de întărire chimică
SK	Chemicky tuhnuce lepidlo
SL	Lepilo za kemično strjevanje
ES	Adhesivo de curado químico
SV	Kemiskt hårdande lim
TR	Kimyasal kür yapıştırıcı
RU	Клей химического отверждения
UA	Клей для хімічного затвердіння

INSTRUCTION FOR USE

DESCRIPTION

dline Chemical Curing Adhesive is polymer-based, bioinert, two component chemical curing adhesive, characterized by excellent adhesion to dentin and enamel. Before use etching of the surface is required.

COMPOSITION

Base: methacrylate mixture 95-100%, initiator 1-5%, inhibitor <1%. Catalyst: methacrylate mixture 95-100%, initiator 1-5%, inhibitor <1%.
dline Chemical Curing Adhesive does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

setting time	<5 minutes
shear bond strength (enamel)	20.7±3.1 Mpa

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Chemical Curing Adhesive helps to restore/improve aesthetic appearance of restorable tooth; helps to restore/maintain dental function of restorable tooth; helps to protect biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For use on dentin and enamel, prior to chemical curing restorative materials to improve retention of a restoration.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

dline Chemical Curing Adhesive should not be used with products containing eugenol as eugenol may disturb polymerization process.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Chemical Curing Adhesive may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth.

INTENDED USER

dline Chemical Curing Adhesive is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental adhesive. There is no need for specific training.

STERILITY

dline Chemical Curing Adhesive is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Chemical Curing Adhesive is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of adhesive is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No consumable components and accessories are supplied with the device.

INSTRUCTION FOR USE

1. Prepare cavity as always. Cleaning with oil-free prophylaxis paste, such as dline Proply Paste is recommended.
2. Apply layer of etch, such as dline Phosphoric Acid Etching Gel on all surface to be etched. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin. After drying, the surface must be chalky and not contaminated prior applying the bonding resin. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
3. Shake bottles before each use.
4. Delicate push on bottle and dispense required amount of adhesive catalyst and base on the mixing pad. Do not force bottle. Mixing ratio by weight: catalyst 1 / base 1. For mixing use plastic spatula and paper pad. Total mixing time is 10 seconds.
5. The working area must be kept dry. Using a fully saturated micro applicator apply layer of adhesive onto etched surface. Total working time is to 1 minute including mixing.
6. After 2 minutes from start of mixing apply the layer of restorative material, such as dline Chemical Curing Composite into cavity and finish restoration.

WARNINGS

Do not use dline Chemical Curing Adhesive for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Chemical Curing Adhesive does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If noticed any changes in the performance of the device, immediately remove applied material with suitable dental instrument. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center.

SHELF-LIFE

Shelf-life of dline Chemical Curing Adhesive is 2 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Chemical Curing Adhesive is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF 14100 3ml adhesive catalyst, 3ml adhesive base

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