
dline Chemical Curing Composite

EN	Chemical curing composite
LT	Cheminiio kietėjimo kompozitas
DE	Chemisch härtendes Komposit
BG	Химически втвърдяващ се композит
CS	Chemicky tuhnutí kompozit
DA	Kemisk hærdende komposit
ET	Keemiliselt kõvenev komposiit
FR	Composite à durcissement chimique
EL	Σύνθετο χημικό πολυμερισμό
HU	Kémiai kötész kompozit
IT	Composito a polimerizzazione chimica
LV	Ķīmiski cietējošs kompozīts
NO	Kjemisk herdende komposit
PL	Kompozyt utwardzany chemicznie
PT	Composto de cura química
RO	Compozit de întărire chimică
SK	Chemicky tuhnutí kompozit
SL	Kemično strjevalni kompozit
ES	Compuesto de curado químico
SV	Kemisk härdningskomposit
TR	Kimyasal sertleşen kompozit
RU	Композит химического отверждения
UA	Хімічний твердіючий композит

INSTRUCTION FOR USE EN

DESCRIPTION

dline Chemical Curing Composite is chemical curing, bioinert restorative material based on methacrylate mixture and radiopaque glass particles.

COMPOSITION

Base: dental glass grinded 70-80%, methacrylate mixture 20-50%, silicon dioxide 1-5%, opacifier <1%, inhibitor <1%, initiator <1%, stabilizer <1%.

Catalyst: dental glass grinded 70-80%, methacrylate mixture 20-50%, silicon dioxide 1-5%, opacifier <1%, inhibitor <1%, initiator <1%, stabilizer <1%.

dline Chemical Curing Composite does not contain medicinal substance, including human blood or plasma derivate; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

curing	chemical
flexural strenght	109.6±8.3 MPa
water sorption	17.20 µg/mm ³
water solubility	4.05 µg/mm ³

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Chemical Curing Composite restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For restoration of class III and V cavities, selected cavities of class IV and for limited use in class I restorations in premolars.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

dline Chemical Curing Composite should not be used with products containing eugenol as eugenol may disturb polymerization process.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Chemical Curing Composite may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body - mouth. Tissues or body fluids contacted by the device - tooth, oral mucosa, saliva.

INTENDED USER

dline Chemical Curing Composite is developed for professional use in dentistry only. Its user only licensed doctor/technician who has knowledge how to use common dental composites. There is no need for specific training.

STERILITY

dline Chemical Curing Composite is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Chemical Curing Composite is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of composite is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as plastic spatulas are supplied with the device.

INSTRUCTION FOR USE

1. Prepare cavity as always. Cleaning with oil-free prophylaxis paste, such as dline Prophyl Paste is recommended. For deep cavities use calcium hydroxide liner or glass ionomer base lining cement, such as dline Glass Ionomer Base Lining Cement.
2. Apply layer of etch, such as dline Phosphoric Acid Etching Gel on all surface to be etched, follow manufacturer's instruction for use. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin.
3. Apply a thin layer of adhesive, such as dline Chemical Curing Adhesive immediately onto etched surface, follow manufacturer's instruction for use.
4. Mix dline Chemical Curing Composite catalyst paste and base paste just before application. Mixing ratio by weight: catalyst 1 / base 1. For mixing use plastic spatula and paper pad. NOTE! Do not use metal spatula for mixing pastes. Total mixing time is 30 seconds.
5. The working area must be kept dry.
6. Apply the mixed material into cavity. Total working time is 2 minutes including mixing.
7. After 6 minutes from start of mixing restoration is ready for polishing and finishing.
8. In order to improve surface and to seal micro cracks and surface imperfections of the finished composite surface use sealant, such as dline Light Curing Pit And Fissure Sealant.

WARNINGS

Do not use dline Chemical Curing Composite for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Chemical Curing Composite does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If noticed any changes in the performance of the device, immediately remove applied material with suitable dental instrument. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center.

SHELF-LIFE

Shelf-life of dline Chemical Curing Composite is 2 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-10°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Chemical Curing Composite is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF 13007 7g catalyst paste, 7g base paste, 3ml adhesive dline Chemical Curing Adhesive catalyst, 3ml adhesive dline Chemical Curing Adhesive base, 3ml etching liquid dline Phosphoric Acid Etching Gel, 5psc. spatulas, 1 mixing pad

REF 13015 15g catalyst paste, 15g base paste, 3ml adhesive dline Chemical Curing Adhesive catalyst, 3ml adhesive dline Chemical Curing Adhesive base, 5ml etching liquid dline Phosphoric Acid Etching Gel, 10psc. spatulas, 2 mixing pads

REF 13030 30g catalyst paste, 30g base paste, 5ml adhesive dline Chemical Curing Adhesive catalyst, 5ml adhesive dline Chemical Curing Adhesive base, 10ml etching liquid dline Phosphoric Acid Etching Gel, 20psc. spatulas, 2 mixing pads

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