# Qline

## Dual Curing Adhesive Resin Cement

EN	Dual curing adhesive resin cement
LT	Dvigubo kietėjimo dervinis surišimo cementas
DE	Dualhärternder adhäsiver Kunststoffzement
BG	Цимент с двойно втвърдяваща се адхезивна смола
CS	Dvojitě tuhnoucí lepicí pryskyřný cement
DA	Dobbelt hærdende klæbende harpikscement
ET	Kahekordsekt kövenev liimvagu tsement
FR	Ciment résine adhésif à double durcissement
EL	Τομέντο ρητίνη με διπλή ακλήρωση
HU	Kettős kötésű ragasztógyanta cement
IT	Cemento resinoso adesivo a doppia polimerizzazione
LV	Divreiz cietējošs līmes sveķu cements
NO	Dobbel herdende limharpiks cement
PL	Klej zywiczny o podwójnym utwardzaniu
PT	Cimento de resina adesiva de cura dupla
RO	Ciment cu rășină adezivă cu întărire duală
SK	Dvojz tvrdnúci lepiaci živicový cement
SL	Cement iz lepilne smole z dvojnim strjevanjem
ES	Cemento de resina adhesiva de curado dual
SV	Dubbelhärdande limharts cement
TR	Cift kürklenen yapışkan reçine çimentosu
RU	Клей-цемент двойного отверждения
UA	Цемент з клейової смоли подвійного затвердіння

### INSTRUCTION FOR USE

EN

#### DESCRIPTION

dline Dual Curing Adhesive Resin Cement is a dual curing, bioinert, radiopaque cement, based on methacrylate resins and silica/barium glass filler particles. dline Dual Curing Adhesive Resin Cement is a high strength, highly resistant device with optimal bonding characteristics. The use of automix tips ensures homogenous mixture which can be placed directly into the restoration.

#### COMPOSITION

Base: dental glass grinded based on silica and barium 50-70%, methacrylate mixture 30-50%, conditioner 1-5%, silicon dioxide 1-5%, initiator <1%, opacifier <1%, accelerator 1%, inhibitor <1%. Catalyst: dental glass grinded based on silica and barium 50-70%, methacrylate mixture 30-50%, silicon dioxide 1-5%, coinitiator <1%, stabilizer <1%, photoinitiator <1%, inhibitor <1%, pigment <1%.

dline Dual Curing Adhesive Resin Cement does not contain medicinal substance, including human blood or plasma derivate; tissues or cells, or their derivates, of human origin; tissues or cells of animal origin, or their derivates, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

#### PERFORMANCE CHARACTERISTICS

dual curing	light and chemical
flexural strength	115.4±8.4 MPa
water sorption	24.28±0.45 µg/mm <sup>3</sup>
water solubility	0.64±0.21 µg/mm <sup>3</sup>
film thickness	8.8 µm

#### INTENDED PURPOSE AND CLINICAL BENEFITS

dline Dual Curing Adhesive Resin Cement restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

#### CLINICAL INDICATIONS

- For the permanent cementation of indirect restorations including ceramic, composite and metal-based inlays, onlays, crowns and bridges;
- For the permanent cementation of posts.

#### CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

#### RESTRICTIONS TO COMBINATIONS

dline Dual Curing Adhesive Resin Cement should not be used with products containing eugenol, hydrogen peroxide as they may disturb polymerization process.

#### UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Dual Curing Adhesive Resin Cement may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract).

#### RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

#### PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

#### INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body - mouth. Tissues or body fluids contacted by the device - tooth.

#### INTENDED USER

dline Dual Curing Adhesive Resin Cement is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental cements. There is no need for specific training.

#### STERILITY

dline Dual Curing Adhesive Resin Cement is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

#### USE ENVIRONMENT

dline Dual Curing Adhesive Resin Cement is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of cement is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

#### CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as automix tips, intraoral-root canal tips are supplied with the device.

#### INSTRUCTION FOR USE

##### PREPARATION:

1. Remove the temporary restoration.
2. Clean the surface, rinse with water and dry with air. Avoid over drying. The surface must remain glossy.
3. Try in the permanent restoration to ensure proper fit.
4. Apply layer of etch, such as dline Phosphoric Acid Etching Gel on all surface to be etched, follow manufacturer's instruction for use. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin. After drying, the surface appears chalky and must not be contaminated prior to applying the adhesive. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
5. Apply a thin layer of light curing adhesive, such as dline Light Curing Adhesive immediately onto etched surface or dual cure adhesive for post cementation, follow manufacturer's instruction for use.
6. Remove dline Dual Curing Adhesive Resin Cement syringe cap.
7. Promptly and carefully attach the automix tip to the syringe.
8. Test flow of material from automix tip before using.
9. Squeeze out a 1-2mm of paste onto a mixing pad to equal the base paste and catalyst paste in the automix tip. Do not use the paste which has been squeezed out.

##### PLACEMENT:

###### INLAYS/ONLAYS:

1. Dispense the cement directly into the restoration. Cover all surfaces.
2. Seat the restoration firmly onto the preparation. Stabilize until cement is fully set.

###### CROWNS/BRIDGES:

1. Dispense the cement directly into the restoration. Cover all surfaces.
2. Seat the restoration firmly onto the preparation. Let the cement flow from all margins. Stabilize until cement is fully set.

###### POSTS:

1. Mount the intraoral-root canal tip onto the automix tip and dispense the cement directly into the canal.
2. Seat the post carefully and rotate slightly to avoid trapping air. Let the cement flow from the canal space. Light cure for 10 seconds to fix the position of the post. Use LED polymerization lamp with light intensity 1200mW/cm<sup>2</sup>. Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use. Stabilize until cement is fully set.

##### FINISHING:

1. After the restoration is properly seated and while the cement is in rubbery state remove all excess cement with scaler or explorer. Rubbery state can be achieved by light cure for approx. 2-3 seconds or by allowing the cement to self-cure for approx. 30 seconds - 1 minute after application.
2. Light cure for 20-30 seconds all surfaces and margins after removal of all excess cement. Use LED polymerization lamp with light intensity 1200mW/cm<sup>2</sup>. Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use. Finish restoration.

##### NOTE!

- a) Working time at 23°C is about 1 minute 30 seconds to 2 minutes. Chemical setting time at 37°C is about 2 minutes to 2 minutes 30 seconds. Light curing time per surface is 20-30 seconds.

##### WARNINGS

Tips are single use only. Do not use in direct contact with pulp tissue. As dline Dual Curing Adhesive Resin Cement highly adhesive, clean all instruments before cement set on. The setting time depends on the ambient and oral temperature. Increased temperature will reduce the working/setting time. Vice versa decreased temperature will prolong the working/setting time. The cement is sensitive to light. The working time significantly reduces during application under intensive lighting. When removing dline Dual Curing Adhesive Resin Cement from refrigerator after storage, leave the material at ambient temperature for some time until warm-up. After the desired amount of material extruded, immediately close the syringe, so that the material is not unlighted. Do not use dline Dual Curing Adhesive Resin Cement for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Dual Curing Adhesive Resin Cement does not emit radiation and does not cause any electromagnetic interferences.

##### PRECAUTIONS!

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

##### SHELF-LIFE

Shelf-life of dline Dual Curing Adhesive Resin Cement is 2 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

##### STORAGE

Keep product tightly closed in dry well-ventilated place at 2-8°C (preferable in refrigerator). Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

##### DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

## VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Dual Curing Adhesive Resin Cement is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

## MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

## VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

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