Glass Ionomer Luting Cement

EN	Glass ionomer luting cement
LT	Stiklo jonomerinis cementas cementavimui
DE	Glasionomer-Befestigungszement
BG	Стъклоно -иономерен цемент за цементиране
CS	Skloionomerní tmelovací cement
DA	Glasionomer -cementering af cement
ET	Klaasionomeerist liimtsement
FR	Ciment de scellement verre ionomère
EL	Τσιμέντο ιοντομερούς γυαλιού
HU	Üvegionomer tömítő cement
IT	Cemento cementizio vetroionomerico
LV	Stikla jonomēru blīvējošais cements
NO	Glasionomer luting sement
PL	Cement glasonomerowy do osadzania
PT	Cimento de cimentação de ionómero de vidro
RO	Ciment ionizant de sticlă ionomer
SK	Skloionomerný cement na tmlenie
SL	Stekloionomerni cement za cementiranje
ES	Cemento de fijación de ionómero de vidrio
SV	Glasjonomer som luter cement
TR	Cam iyonomer yapıştırma simanı
RU	Стеклоиономерный цемент для фиксации
UA	Склоіономерний цемент для цементування

EN ISO 9917-1 4.2a Glass Polyalkenoate Cement

INSTRUCTION FOR USE EN

DESCRIPTION

dline Glass Ionomer Luting Cement is a chemical curing, bioinert, glass ionomer fixing cement with high adhesion to dentin, enamel and non-precious alloys. Has low film thickness, good translucency and excellent biocompatibility.

COMPOSITION

Fluoroaluminosilicate glass 75-95%, water 20-25%, polyacrylic acid 5-25%, excipient 0-2%, pigment 0-1%.

dline Glass Ionomer Luting Cement does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

compressive strength	126.1 MPa
film thickness	15.4 µm
lead content	0.17 mg/kg

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Glass Ionomer Luting Cement restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For cementing metal and ceramic crowns, bridges, inlays, onlays.
- For cementing orthodontic brackets and bands.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Glass Ionomer Luting Cement may cause allergic or irritation reactions (eye).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OR BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth, oral mucosa, saliva.

INTENDED USER

dline Glass Ionomer Luting Cement is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental cements. There is no need for specific training.

STERILITY

dline Glass Ionomer Luting Cement is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Glass Ionomer Luting Cement is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of cement is suitable for single

use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as measuring scoop, mixing pads, are supplied with the device.

INSTRUCTION FOR USE

1. Prepare surface as always. For deep cavities use calcium hydroxide liner.
2. Mixing ratio: 1 level spoon of powder / 2 drops of liquid. Lightly tap the bottle of powder against the hand before dispensing. Do not shake or invert. Hold liquid bottle vertically upright and squeeze gently when dispense a liquid.
3. For mixing use metal or plastic spatula and glass slab or paper mixing pad. Divide powder into two parts. Mix the first half of powder with liquid for 10 seconds, and then add remainder. Mix until exact consistency is reached. Total mixing time is 30 seconds. Do not add powder in small portions.
4. The working area must be kept dry. Apply the mixed cement with suitable instrument to both restoration and prepared tooth immediately after mixing. Total working time is 2 minutes including mixing. Excess cement may be removed immediately or after cement has set. Do not attempt to remove excess cement in the phase of setting (1.5-6 minutes from start of mix) as adhesion may be impaired.
5. After 8 minutes from start of mixing restoration is ready for finishing.

WARNINGS

As dline Glass Ionomer Luting Cement highly adhesive, clean all instruments before cement set on. Higher amount of powder or higher temperature shorten the working/setting time and lower temperature, or higher amount of liquid prolong the working/setting time. Do not use dline Glass Ionomer Luting Cement for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Glass Ionomer Luting Cement does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Wash hands thoroughly after handling. Avoid contact with eyes. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If during the use of the product, after all the components are well mixed, product consistency is non-homogenous, or product does not harden or does not adhere to the applied surface while setting or any by-products/phases are released during the time of curing, or sudden acute pain occur on application site, or if any other abnormal behavior of the product is noticed while manipulating the device, that is not mentioned above, discontinue to use immediately. Remove the restoration from the tooth cavity with suitable dental instrument do not let the product to be swallowed. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center. Collect all available remaining supplies, do not use them again and keep them out of reach in a safe place until further notice. Contact the manufacturer immediately and report of any noticed changes in the performance of the product.

SHELF-LIFE

Shelf-life of dline Glass Ionomer Luting Cement is 3 years of powder / 5 years of liquid from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 5-25°C. Protect from moisture, direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Glass Ionomer Luting Cement is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance is available through manufacturer's website www.dline.ee/sscp/ until European Database on Medical Devices (EUDAMED) comes online.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF 16400 20g powder, 13g liquid, measuring scoop, mixing pad

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