dline Light Curing Adhesive

EN	Light curing adhesive
LT	Šviesoje kietėjantis surišėjas
DE	Lichthärtender Klebstoff
BG	Лекo втвърдяващо лепило
CS	Světlem tuhnutí lepidlo
DA	Let hærdende klæbemiddel
ET	Kergelt kõvenev liim
FR	Adhésif photopolymérisable
EL	Κόλλη ελαφρώς ακλήρυνσης
HU	Fényre kötő ragasztó
IT	Adesivo fotopolimerizzabile
LV	Gaismas cietēšanas lieme
NO	Letthærdende lim
PL	Klej światłoutwardzalny
PT	Adesivo fotopolimerizável
RO	Adeziv de întărire ușoară
SK	Svetlo tuhnuce lepidlo
SL	Lepilo za utrjevanje svetlobe
ES	Adhesivo fotopolimerizable
SV	Lätthårdande lim
TR	Işıkla sertleşen yapıştırıcı
RU	Светоотверждаемый клей
UA	Світлотвердуючий клей

INSTRUCTION FOR USE

DESCRIPTION

dline Light Curing Adhesive is a 5th generation light curing, bioinert, single component one-step adhesive, suitable for use with wet bonding technique. Etching of the surfaces is required.

COMPOSITION

Methacrylate mixture 55-95%, acetone 10-30%, conditioner 1-10%, coinitiator <5%, accelerator <1%, inhibitor <1%, photoinitiator <1%.

dline Light Curing Adhesive does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

light curing	20 s
shear bond strength (enamel)	21.8±2.4 MPa

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Light Curing Adhesive helps to restore/improve aesthetic appearance of restorable tooth; helps to restore/maintain dental function of restorable tooth; helps to protect biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- **Used to bond dentin and enamel to polymer-based restorative materials, comomers, ceramic surfaces to improve retention of a restoration.**

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

dline Light Curing Adhesive should not be used with products containing eugenol as eugenol may disturb polymerization process.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Light Curing Adhesive may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth.

INTENDED USER

dline Light Curing Adhesive is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental adhesive. There is no need for specific training.

STERILITY

dline Light Curing Adhesive is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Light Curing Adhesive is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of adhesive is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as micro applicators, cup for application, are supplied with the device.

INSTRUCTION FOR USE

1. Prepare cavity as always. Cleaning with oil-free prophylaxis paste, such as dline Proply Paste is recommended. For deep cavities use calcium hydroxide liner or glass ionomer based cement.
2. Apply layer of etch, such as dline Phosphoric Acid Etching Gel on all surface to be etched. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin. After drying, the surface must be chalky and do not contaminated prior to applying the bonding resin. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
3. Shake bottle before each use.
4. Delicate push on bottle and dispense small amount of adhesive into cup. Do not force bottle.
5. Using a fully saturated micro applicator apply 2-3 layers of adhesive onto etched surface.
6. Dry with air for 10-20 seconds to dry all solvent from surface. If surface is not glossy, apply additional layers of adhesive and dry.
7. Curing with light for 20 seconds. Use LED polymerization lamp with light intensity 1200mW/cm². Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use.
8. Apply the layer of light curing restorative material, such as dline Light Curing Zirconium Reinforced Composite, dline Light Curing Nano Hybrid Composite or dline Light Curing Nano Flowable Composite into cavity and finish restoration.

WARNINGS

After the desired amount of adhesive extruded, immediately close the bottle. dline Light Curing Adhesive is highly flammable liquid. Do not use dline Light Curing Adhesive for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Light Curing Adhesive does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If noticed any discoloration of the patient's tooth tissue during use, strong pungent odor, uneven/non-homogenous consistency, or patient suffers from sudden acute pain on application site, or if any other abnormal behavior of the product noticed while manipulating the device, that was not mentioned above, discontinue to use immediately. Rinse applied product with plenty of water, use suction system and do not let the product enter the mouth and/or to be swallowed. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center. Collect all available remaining supplies, do not use them again and keep them out of reach in a safe place until further notice. Contact the manufacturer immediately and report of any noticed changes in the performance of the product.

SHELF-LIFE

Shelf-life of dline Light Curing Adhesive is 3 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Light Curing Adhesive is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance is available through manufacturer's website www.dline.ee/sscp/ until European Database on Medical Devices (EUDAMED) comes online.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF 14000

5ml bottle, 5pcs micro applicators, cup

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