dline Light Curing Compomer Liner

EN	Light curing compomer liner
LT	Sviesoje kietėjantis kompomero pamušalas
DE	Lichthärtender Compomer-Liner
BG	Леко втвърдяваща се компомерна подложка
CS	Světlem tuhnoucí kompomerní vložka
DA	Let hærdende kompomerfor
ET	Valguskõvastuv kompomereeri vooder
FR	Liner compomère photopolymérisable
EL	Ελαφριά επένδυση κομπόμερ
HU	Fényre kötő kompomér béles
IT	Liner in compomero fotopolimerizzabile
LV	Gaismas cietēšanas kompomēre
NO	Letthærdende kompomerfor
PL	Kompomer światłoutwardzalny
PT	Forro de compômero fotopolimerizável
RO	Căptuşeală pentru întărire uşoară
SK	Svetlom tuhnuća vložka z kompoméru
SL	Komponentna podloga za svetlo utrjevanje
ES	Revestimiento de compómtero fotopolimerizable
SV	Lätthårdande kompomerfoder
TR	Işıkla sertleşen kompomér astar
RU	Светоотверждаемый компомерный лайнер
UA	Компомерний вкладиш із легким затвердінням

INSTRUCTION FOR USE

EN

DESCRIPTION

dline Light Curing Compomer Liner is light curing, bioinert, radiopaque compomer liner based on polyacid-modified composite resins and glass ionomer filler. Such composition allows compomer expansion partially compensate curing contraction of composite materials. Thanks to the glass ionomer content compomer doubles the adhesion to dentin. dline Light Curing Compomer Liner prevents crack formations, minimizes micro fractures and secondary caries. Material could be used both as a usual liner (coating the floor of the cavity) and to coat the entire cavity (stress-breaking function).

COMPOSITION

Dental glass reactive 50-70%, methacrylate mixture 20-30%, modifying agent 1-5%, photoinitiator <1%, cointiator <1%, opacifier <1%, stabilizer <1%, inhibitor <1%.

dline Light Curing Compomer Liner does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

light curing	20s
depth of cure	2.51±0.11 mm
flexural strength	34.3±2.8 MPa

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Light Curing Compomer Liner helps to restore/improve aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For lining under polymer-based restorative materials, amalgams and porcelain restorations.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

dline Light Curing Compomer Liner should not be used with products containing eugenol as eugenol may disturb polymerization process.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Light Curing Compomer Liner may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth.

INTENDED USER

dline Light Curing Compomer Liner is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental liners. There is no need for specific training.

STERILITY

dline Light Curing Compomer Liner is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or

calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Light Curing Compomer Liner is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of liner is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as application tips, are supplied with the device.

INSTRUCTION FOR USE

1. Prepare cavity as always. Cleaning with oil-free prophylaxis paste, such as dline Prophyl Paste is recommended. Remove all cleaning residues with water spray and dry with air. Avoid contamination with saliva during applications.
2. In deep cavities use calcium hydroxide liner.
3. Apply layer of etch, such as dline Phosphoric Acid Etching Gel to all surface to be etched, follow manufacturer's instruction for use. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse surface with water and dry with air. Avoid over drying dentin.
4. Apply a layer of adhesive, such as dline Light Curing Adhesive immediately onto etched surface, follow manufacturer's instruction for use. The tooth cavity must be kept dry.
5. dline Light Curing Compomer Liner could be applied without adhesive. In this case, adhesion is similar to that obtained with usual cements (glass ionomer cement, phosphate cement).
6. Using application tip apply layer of dline Light Curing Compomer Liner in a thickness of 1mm.
7. Cure with light for 20 seconds at least. Use LED polymerization lamp with light intensity 1200mW/cm². Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use.
8. Apply the layer of composite, such as dline Light Curing Zirconium Reinforced Composite, dline Light Curing Nano Hybrid Composite, dline Light Curing Nano Flowable Composite, follow manufacturer's instruction for use, and finish restoration.

WARNINGS

After the desired amount of material extruded, immediately remove application tip and close the syringe cap, so that the material is not unlighted. The material is sensitive to light. Avoid too long manipulation time under intensive lighting. Compomer adheres to metal instruments and should be removed with a paper towel before it has set. Place the lights as close as possible to the material. Use of a matrix is recommended. Do not use dline Light Curing Compomer Liner for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Light Curing Compomer Liner does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If during the use of the product noticed any abnormal product performance characteristics: non-homogenous, non-flowable, uneven consistency, product does not harden or does not adhere on the applied surface while curing with light or any by-products/phases are released during curing, or sudden acute pain occur on application site, or if any other abnormal behavior of the product noticed while manipulating the device, that is not mentioned above, discontinue to use immediately. Remove the restoration from the tooth cavity with suitable dental instrument do not let the product to be swallowed. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center. Collect all available remaining supplies, do not use them again and keep them out of reach in a safe place until further notice. Contact the manufacturer immediately and report of any noticed changes in the performance of the product.

SHELF-LIFE

Shelf-life of dline Light Curing Compomer Liner is 4 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Light Curing Compomer Liner is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance is available through manufacturer's website www.dline.ee/sscp/ until European Database on Medical Devices (EUDAMED) comes online.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

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