dline Light Curing Nano Flowable Composite

EN	Light curing nano flowable composite
LT	Šviesoje kietėjantis nano takus kompozitas
DE	Lichthärtendes nano-fließfähiges Komposit
BG	Леко втвърдяващ се нанотечен композит
CS	Světlem tuhnoucí nano tekutý kompozit
DA	Let hærdende nano flydende komposit
ET	Valguskõvastuv nano voolav komposiit
FR	Composite nano-fluide photopolymérisable
EL	Ελαφρύ σκληρυνση νανο -ρευστό σύνθετο
HU	Fényben kötő nano folyékony kompozit
IT	Composito fluido fotopolimerizzabile nano
LV	Gaisā cietējošs nano plūstošs kompozīts
NO	Lysherdende nanoflytbar komposit
PL	Światłoutwardzalny nano płynny kompozyt
PT	Composto nanofluido fotopolimerizável
RO	Compozit nano fluid cu fotopolimerizare
SK	Svetlom tuhnutí nano tekutý kompozit
SL	Nano tekoči kompozit, polimeriziran s svetlobo
ES	Composite nano fluido fotopolimerizable
SV	Lätthårdande nanoflytbar komposit
TR	İşikla sertleşen nano akışkan kompozit
RU	Светоотверждаемый нано-текучий композит
UA	Світлотвердіючий нано текучий композит

INSTRUCTION FOR USE EN

DESCRIPTION

dline Light Curing Nano Flowable Composite is a light curing nano flowable, bioinert, radiopaque composite under the Vita* shades. dline Light Curing Nano Flowable Composite is a high aesthetic, highly resistant, superior polishability product with optimal flow characteristics.

COMPOSITION

Dental glass grinded 50-70%, methacrylate mixture 30-40%, silicon dioxide 1-5%, coinitiator <1%, photoinitiator <1%, stabilizer <1%, inhibitor <1%, opacifier <1%, pigment 1%.

dline Light Curing Nano Flowable Composite does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

light curing	20-30s
depth of cure	3.24±0.03 mm
flexural strength	110.2±4.1 MPa
water sorption	17.10±0.20 µg/mm ³
water solubility	0.00±0.00 µg/mm ³

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Light Curing Nano Flowable Composite restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For restorations of class III, IV and V cavities; root surface caries restorations;
- For sealing pits and fissures;
- For initial placement in class I and II cavities.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

dline Light Curing Nano Flowable Composite should not be used with products containing eugenol as eugenol may disturb polymerization process.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Light Curing Nano Flowable Composite may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth, saliva.

INTENDED USER

dline Light Curing Nano Flowable Composite is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental composites. There is no need for specific training.

STERILITY

dline Light Curing Nano Flowable Composite is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular

maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Light Curing Nano Flowable Composite is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of composite is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as application tips, are supplied with the device.

INSTRUCTION FOR USE

CAVITY PREPARATION:

1. Prepare cavity as always.
2. Clean the surface with oil-free prophylaxis paste, such as dline Prophy Paste.
3. For deep cavities use calcium hydroxide liner or glass ionomer base lining cement, such as dline Glass Ionomer Base Lining Cement.

ETCHING, BONDING:

1. Apply layer of etch, such as dline Phosphoric Acid Etching Gel to surface to be etched. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin.
2. Apply a layer of adhesive, such as dline Light Curing Adhesive immediately onto etched surface, follow manufacturer's instruction for use.

SYRINGE PREPARATION:

1. Remove syringe cap.
2. Promptly and carefully attach the dispensing tip to the syringe.
3. Test flow of materials from tip before using intraorally.

PLACEMENT OF dline Light Curing Nano Flowable Composite:

1. Before bringing the syringe to the mouth, remove the air from the dispensing tip. To remove air from the tip, with the tip pointing upwards, gently push forward the syringe plunger. If the air is still inside the dispensing tip, air bubbles may be removed at the time of injection.
2. Delicate push on plunger and apply layer of material into the cavity. Do not force plunger.
3. Do not apply layers more than 2 mm deep.
4. Light cure for 20-30 seconds (depends on layer deep). Use LED polymerization lamp with light intensity 1200mW/cm² in full mode, not ramp or pulse mode. Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use. Finish restoration.

WARNINGS

After the desired amount of material extruded, immediately remove application tip and close the syringe cap, so that the material is not unlighted. The material is sensitive to light. Avoid too long manipulation time under intensive lighting. Do not use dline Light Curing Nano Flowable Composite for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Light Curing Nano Flowable Composite does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

SHELF-LIFE

Shelf-life of dline Light Curing Nano Flowable Composite is 4 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Light Curing Nano Flowable Composite is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

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* Registered trademark of the Vita Zahnfabrik H.Rauter GmbH & Co. KG, Bad Sackingen, Germany.