Phosphoric Acid Etching Gel

EN	Phosphoric acid etching gel
LT	Fosforo rūgšties išdinimo gelis
DE	Ätzelgel mit Phosphorsäure
BG	Гел за ецване на фосфорна киселина
CS	Gél pro leptání kyselinou fosforečnou
DA	Phosphorsyre ætsende gel
ET	Fosforhappe söövitat geel
FR	Gel de mordantçage à l'acide phosphorique
EL	Τζελ χύραξης φωσφορικού οξέος
HU	Foszforsav marató gél
IT	Gel mordenzante con acido fosforico
LV	Fosforskābes kodināšanas želeja
NO	Fosforsyre etsegel
PL	Żel do wytrawiania kwasem fosforowym
PT	Gel de ácido fosfórico
RO	Gel de gravat cu acid fosforic
SK	Gél na leptanie kyselinou fosforečnou
SL	Gel za jedkanje fosforne kisline
ES	Gel de grabado de ácido fosfórico
SV	Fosforsyra etsgel
TR	Fosforik asit aşındırma jeli
RU	Гель для травления с фосфорной кислотой
UA	Гель для травлення фосфорної кислоти

INSTRUCTION FOR USE EN

DESCRIPTION

dline Phosphoric Acid Etching Gel is 37% thixotropic phosphoric acid gel with optimum level of viscosity. Syringe presentation with special tips allows placement of the etching gel in any desired location.

COMPOSITION

Diluent 40-60%, phosphoric acid 37%, silicon dioxide 10-20%, thickener 5-10%, pigment <1%.

dline Phosphoric Acid Etching Gel does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

pH	<1
----	----

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Phosphoric Acid Etching Gel is intended to be used for tooth surface preparation procedures.

CLINICAL INDICATIONS

- Tooth enamel and / or dentin etching.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Phosphoric Acid Etching Gel may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth.

INTENDED USER

dline Phosphoric Acid Etching Gel is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental etchants. There is no need for specific training.

STERILITY

dline Phosphoric Acid Etching Gel is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Phosphoric Acid Etching Gel is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of etchant is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as application tips, are supplied with the device.

INSTRUCTION FOR USE

PREPARATION:

1. Prepare surface as always. The working area must be kept dry.

2. If the surface to be etched is not freshly exposed, cleaning with oil-free prophylaxis paste such as dline Propyl Paste is recommended.
3. Attach application tip to the syringe.
4. Flow small amount of gel on mixing pad.

ETCHING:

1. Delicate push on plunger and apply layer of etch on all surface to be etched. Do not force plunger.
2. Leave etch gel in place for 15 seconds (dentine), 30 seconds (enamel).
3. Rinse with water and dry with air.
4. Verify etched surface is clean and dry.
5. Use an adhesive such as dline Light Curing Adhesive.

WARNINGS

After the desired amount of material extruded, immediately remove application tip and close the syringe cap. Do not use dline Phosphoric Acid Etching Gel for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Phosphoric Acid Etching Gel does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. Subsequently, have the casualty drink 1 - 2 glasses of water immediately. Do not try to neutralize with bases, no charcoal. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

Precautions to be taken in the event of changes in the performance of the device:

If noticed any discoloration and/or increased gloss of the patient's tooth tissue during use, sudden acute pain, strong pungent odor, foaming, uneven/run-off application or any other abnormal behaviour of the product, while manipulating the device, that is not mentioned above, discontinue to use immediately. Rinse applied gel with plenty of water, use suction system and do not let the product enter the mouth and/or to be swallowed. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center. Continue to collect all available remaining supplies, do not use them again and keep them out of reach in a safe place until further notice. Contact the manufacturer immediately and report of any noticed changes in the performance of the product.

SHELF-LIFE

Shelf-life of dline Phosphoric Acid Etching Gel is 3 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Phosphoric Acid Etching Gel is safe and performs as intended if it is used in accordance to manufacturer's instruction for use.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

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