Silver-Glass Ionomer Cement

EN	Silver-glass ionomer cement
LT	Sidabro-stiklo jonomerinis cementas
DE	Silber-Glas-Ionomerzement
BG	Сребърно-стъклен йонимерен цемент
CS	Ionomerní cement ze stříbrného skla
DA	Sølvglas ionomer cement
ET	Hõbedast klaasist ionomeertsement
FR	Ciment verre-argent ionomère
EL	Τσιμέντο ιοντομερούς από ασήμι γυαλί
HU	Ezüst-üveg ionomer cement
IT	Cemento vetroionomerico argento
LV	Sudraba stikla jonomēra cements
NO	Sølvglass ionomer sement
PL	Cement srebrno-glasjonomerowy
PT	Cimento de ionómero de prata-vidro
RO	Ciment ionomer de sticlă argintie
SK	Ionomerný cement striebro-sklo
SL	Srebro-stekloionomerni cement
ES	Cemento de ionómero de vidrio y plata
SV	Jonomercement av silverglas
TR	Gümüş-sam iyonomer siman
RU	Цемент стеклоиономерный серебрястый
UA	Іономерний цемент зі срібля

EN ISO 9917-1 4.2b & 4.2c Glass Polyalkenoate Cement

INSTRUCTION FOR USE EN

DESCRIPTION

dline Silver-Glass Ionomer Cement is a chemical curing, bioinert, fast-setting, radiopaque, non-translucent silver reinforced glass ionomer cement with excellent thermal expansion and high strength and hardness.

COMPOSITION

Fluoroaluminosilicate glass 60-80%, silver alloy <25%, zirconium oxide <25%, polyacrylic acid <15%, water <15%, opacifier <10%, excipient <5%, pigment <1%. dline Silver-Glass Ionomer Cement does not contain medicinal substance, including human blood or plasma derivate; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

compressive strength	232.0 MPa
----------------------	-----------

INTENDED PURPOSE AND CLINICAL BENEFITS

dline Silver-Glass Ionomer Cement restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For restoration of class I and II cavities; base under posterior composite restorations and amalgams; core build-ups; crowns and amalgams repairs; deciduous teeth fillings.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, dline Silver-Glass Ionomer Cement may cause allergic or irritation reactions (eye).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OR BODY FLUIDS

Part of the body - mouth. Tissues or body fluids contacted by the device - tooth, oral mucosa, saliva.

INTENDED USER

dline Silver-Glass Ionomer Cement is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental cements. There is no need for specific training.

STERILITY

dline Silver-Glass Ionomer Cement is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

dline Silver-Glass Ionomer Cement is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of cement is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as measuring scoop, mixing pads, are supplied with the device.

INSTRUCTION FOR USE

1. Prepare surface as always. For deep cavities use calcium hydroxide liner. For the best chemical bonding to the teeth, apply conditioner, such as dline Conditioner For Glass Ionomer Cements on prepared surface for 10 sec. Wash with water and dry with air.
2. Mixing ratio: 2 level spoons of powder / 1 drop of liquid. Lightly tap the bottle of powder against the hand before dispensing. Do not shake or invert. Hold liquid bottle vertically upright and squeeze gently when dispense a liquid.
3. For mixing use metal or plastic spatula and glass slab or paper mixing pad. Divide powder into two parts. Mix the first half of powder with liquid for 10-15 seconds, and then add remainder. Mix until exact consistency is reached. Total mixing time is 30 seconds. Do not add powder in small portions.
4. The working area must be kept dry. Apply the mixed cement with suitable instrument into cavity. Total working time is 2 minute including mixing.
5. After 8 minutes from start of mixing restoration is ready for finishing.

WARNINGS

As dline Silver-Glass Ionomer Cement highly adhesive, clean all instruments before cement set on. Higher amount of powder or higher temperature shorten the working/setting time and lower temperature, or higher amount of liquid prolong the working/setting time. Do not use dline Silver-Glass Ionomer Cement for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. dline Silver-Glass Ionomer Cement does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Wash hands thoroughly after handling. Avoid contact with eyes. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

SHELF-LIFE

Shelf-life of dline Silver-Glass Ionomer Cement is 5 years of powder / 5 years of liquid from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 5-25°C. Protect from moisture, direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

dline Silver-Glass Ionomer Cement is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF 16200 10g powder, 3g liquid, measuring scoop, mixing pad

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